

DNA concentration by precipitation in ethanol

based on protocol found at: <https://www.auburn.edu/~santosr/protocols/NaAcEtOHDNAPrecip.pdf>

1. Adjust the salt concentration, for example, with sodium acetate (0.3 M, pH 5.2, final concentration) or ammonium acetate (2.0–2.5 M, final concentration).
2. Add 0.6–0.7 volumes of room-temperature isopropanol to the DNA solution and mix well.
3. Centrifuge the sample immediately at 10,000–15,000 x g for 15–30 min at 4°C
4. Carefully decant the supernatant without disturbing the pellet.
5. Wash the DNA pellet by adding 1–10 ml (depending on the size of the preparation) of room-temperature 70% ethanol. This removes co-precipitated salt and replaces the isopropanol with the more volatile ethanol, making the DNA easier to re-dissolve.
6. Centrifuge at 10,000–15,000 x g for 5–15 min at 4°C
7. Repeat the washing and centrifugation steps once more.
8. Carefully decant the supernatant without disturbing the pellet.
9. Air-dry the pellet for 5–20 min (depending on the size of the pellet).
10. Re-dissolve the DNA in a suitable buffer.