

## Coverslip Immunofluorescence

1. Remove coverslips stored in methanol, and place in a sterile well at a slanted angle to allow the methanol to evaporate (~10-15 min). [As an alternative to methanol sterilization, 70% ethanol can be used for a minimum of 30 min.]
2. After methanol is evaporated, gently tap plate to slide coverslip flat onto the bottom of the well. Desired number of cells can now be plated into the 6/12 well plate as usual.
3. After cells have grown to desired confluence, starve and treat with drugs/growth factors, or just move onto the next step.
4. Wash 2x in PBS.
5. Fix with 4% paraformaldehyde in PBS for 20 min in fume hood. Caution: Formaldehyde is very toxic. Handle with care and make sure all work is done in proper hood.
6. You may either store overnight at 4°C and resume next day after this step, or continue. Wash 2x with PBS, and then permeabilize cells with 0.25% Triton-X 100 in PBS for 5 min.  
*Note: May need to try different detergents for different amounts of times to optimize for your cells + antibodies. Don't go too long on this step however, or cells will lyse.*
7. Wash 2x in PBS.
8. Prepare humidified chamber from any box that can be made airtight.
  - a. Place filter paper on bottom, and add H<sub>2</sub>O to wet filter paper. Pour out excess water.
  - b. Add parafilm on top of filter paper, then pat down.
9. Prepare primary antibody at desired dilution (e.g., 1:500 for Shp2, 1:250 for EGFR) in Odyssey Blocking Buffer. Abs can also be diluted in 3% BSA in PBS, but OBB tends to provide superior results.  
*Note: Centrifuge antibody for 1min at max RPMs in order to settle any particulates.*
10. Use tweezers/forceps to transfer coverslips to the humidified chamber, cells facing up.
  - a. Dot slide on a Kimwipe first to dry before putting into the chamber.
11. Add 100-150 $\mu$ L (125  $\mu$ L for a size 18 coverslip) of antibody solution to each coverslip, just need enough to equally coat the entire surface.
12. Incubate at 37°C for 3 hr. Can either use a water bath or incubator, just something that will not be disturbed. Alternatively, can incubate overnight in the cold room.
13. Transfer coverslips back to 6/12 well plate, and wash 5x in PBS-Tween for approx. 3 min on a rotator/rocker.
14. Prepare secondary antibody at desired dilution (1:750 of antibody that is already diluted 1:2 in glycerol) in OBB. Can also add DNA stain if desired (e.g., 1:2000 of Hoescht).  
*Note: Make sure you use Alexa Fluor or other visible spectrum secondary antibodies, not Licor abs!*  
*Note: Centrifuge antibody for 10 min at max RPMs to settle any particulates.*
15. Transfer coverslips to humidified chamber, dotting on a Kimwipe beforehand. Add 100-150  $\mu$ L of antibody solution to each coverslip to coat the surface.
16. Incubate at 37°C for 1 hr.
17. Wash in 0.1% PBS-Tween20 5x for approx. 3 min on a rotator/rocker.
18. Place a very small drop of mounting media (Prolong Gold Antifade) on a labeled slide, and gently put coverslip onto drop with cells facing down. Allow to dry overnight at room temp in a dark area on a level surface.  
*Note: Mounting media is very viscous. Cut tip off partially on a 200  $\mu$ L tip, and draw up 125  $\mu$ L into cut pipet tip for a total of 6 coverslips. Avoid making bubbles.*
19. After drying overnight, carefully paint a layer of clear nail polish around the edge of the coverslip to prevent mounting media from drying out quickly.