

Transforming competent cells

Before you begin: Be sure the water bath is turned on and set to 42°C.

1. On ice: thaw competent *E. coli* (kept in the -80°C freezer). If using a new 0.5 mL tube of *E. coli* (applies to DH5 α or XL1 Blue), aliquot into 50 μ L aliquots in 1.5 mL Eppendorf tubes. For each transformation, use 50 μ L of DH5 α , XL1 Blue, or Stbl3 cells. Add 1 - 5 μ L DNA to the mixture. [Typically, this means somewhere in the range of 1-100 ng, depending on the source of the DNA, e.g., low yield extraction from filter paper or high concentration stock from a miniprep or maxiprep.] Be careful pipetting and mixing cells as they are fragile; do not pipet up and down.
2. Incubate the mixture of bacteria and DNA on ice with occasional swirling for 30 min.
3. Heat shock bacteria by floating tubes in a 42°C water bath for 40-45 s.
4. Incubate tubes on ice for 2 min.
5. Add 750 μ L SOC medium per sample, place the Eppendorf tube inside a plastic 15 mL culture tube, and incubate for 1 hr with shaking (220 rpm) at 37°C. [If necessary, a 30 min minimum incubation is fine.] During the 1 hr incubation, put ampicillin LB agar plates into the 37°C dry air incubator to warm.
6. Spread 50-300 μ L of bacteria onto a plate using a few glass beads (autoclaved) or a sterile pipet tip. [If needed to plate more cells, you can spread a larger volume of bacteria on the plate. To minimize the total volume of fluid added to the plate, you can spin down the bacteria and re-constitute in a lower volume.]
7. Incubate plates overnight at 37°C. Pick colonies about 12-16 hr later. [Plates with colonies can be kept in the cold room or in a refrigerator dedicated to microbiological work.]